

RN Pipette Tips



Compatibility: Rainin LTS Series

Features

- Imported polypropylene (PP), USP Class VI standards
- Smooth inner wall, reduce liquid residue, ensure the accuracy of liquid suction
- DNA-free, DNase-free, RNase-free, Pyrogen/Endotoxin free
- Ultra-grinding & molding technology, release agent free to ensure better product quality
- Tolerance temperature ranging from -80 °C to 121 °C, stable performance under high temperature and high pressure
- Extra thin-wall design increase product flexibility and improve heat transfer efficiency



RN 20 µL Pipette Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600101	20 µL	Bulk	Clear	1000 pcs/pack, 10 packs/case
8600103	20 µL	Rack packed, Sterile	Clear	96 pcs/rack, 50 racks/case
8600104	20 µL	Low retention, Bulk	Clear	1000 pcs/pack, 10 packs/case
8600106	20 µL	Low retention, Rack packed, Sterile	Clear	96 pcs/rack, 50 racks/case

RN 20 µL Aerosol Filter Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600107	20 µL	Bulk, Filtered	Clear	1000 pcs/pack, 10 packs/case
8600108	20 µL	Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 50 racks/case
8600109	20 µL	Low retention, Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 50 racks/case
8600119	20 µL	Low retention, Bulk, Filtered	Clear	1000 pcs/pack, 10 packs/case

RN 200 µL Pipette Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600201	200 µL	Bulk	Clear	1000 pcs/pack, 20 packs/case
8600203	200 µL	Rack packed, Sterile	Clear	96 pcs/rack, 50 racks/case
8600204	200 µL	Low retention, Bulk	Clear	1000 pcs/pack, 20 packs/case
8600206	200 µL	Low retention, Rack packed, Sterile	Clear	96 pcs/rack, 50 racks/case

RN 200 µL Aerosol Filter Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600207	200 µL	Bulk, Filtered	Clear	1000 pcs/pack, 20 packs/case
8600208	200 µL	Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 50 racks/case
8600209	200 µL	Low retention, Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 50 racks/case
8600219	200 µL	Low retention, Bulk, Filtered	Clear	1000 pcs/pack, 20 packs/case

RN 300 µL Pipette Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600301	300 µL	Bulk	Clear	1000 pcs/pack, 20 packs/case
8600303	300 µL	Rack packed, Sterile	Clear	96 pcs/rack, 40 racks/case
8600304	300 µL	Low retention, Bulk	Clear	1000 pcs/pack, 20 packs/case
8600306	300 µL	Low retention, Rack packed, Sterile	Clear	96 pcs/rack, 40 racks/case

RN 300 µL Aerosol Filter Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600307	300 µL	Bulk, Filtered	Clear	1000 pcs/pack, 20 packs/case
8600308	300 µL	Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 40 racks/case
8600309	300 µL	Low retention, Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 40 racks/case
8600319	300 µL	Low retention, Bulk, Filtered	Clear	1000 pcs/pack, 20 packs/case



RN 1000 µL Pipette Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600401	1000 µL	Bulk	Clear	1000 pcs/pack, 5 packs/case
8600403	1000 µL	Rack packed, Sterile	Clear	96 pcs/rack, 40 racks/case
8600404	1000 µL	Low retention, Bulk	Clear	1000 pcs/pack, 5 packs/case
8600406	1000 µL	Low retention, Rack packed, Sterile	Clear	96 pcs/rack, 40 racks/case

RN 1000 µL Aerosol Filter Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600407	1000 µL	Bulk, Filtered	Clear	1000 pcs/pack, 5 packs/case
8600408	1000 µL	Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 40 racks/case
8600409	1000 µL	Low retention, Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 40 racks/case
8600419	1000 µL	Low retention, Bulk, Filtered	Clear	1000 pcs/pack, 5 packs/case

RN 1200 µL Pipette Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600501	1200 µL	Bulk	Clear	1000 pcs/pack, 5 packs/case
8600503	1200 µL	Rack packed, Sterile	Clear	96 pcs/rack, 40 racks/case
8600504	1200 µL	Low retention, Bulk	Clear	1000 pcs/pack, 5 packs/case
8600506	1200 µL	Low retention, Rack packed, Sterile	Clear	96 pcs/rack, 40 racks/case

RN 1200 µL Aerosol Filter Tips

P/N	Max Capacity	Tip Format	Tips Color	Configuration
8600507	1200 µL	Bulk, Filtered	Clear	1000 pcs/pack, 5 packs/case
8600508	1200 µL	Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 40 racks/case
8600509	1200 µL	Low retention, Rack packed, Filtered, Sterile	Clear	96 pcs/rack, 40 racks/case
8600519	1200 µL	Low retention, Bulk, Filtered	Clear	1000 pcs/pack, 5 packs/case

