

# Screw Cap Micro Tube



- Screw cap tubes are available in a variety of colors with a conical or self-standing bottom
  - Colored-cap versions are also available
  - Manufactured from high-clarity, ultra-pure polypropylene
  - Suitable for a wide range of temperatures-from -80 °C to +120 °C.
- Note: not suitable for immersion in liquid nitrogen
- Ideal for boiling applications
  - Skirted or conical bases
  - Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free



### Screw Cap Micro Tubes Type S

P/N	Capacity	Tube Color	Cap Color	Configuration
801205	0.5 mL	Clear	Clear	100 tubes/pack, 10 packs/case
801207	1.5 mL	Clear	Clear	100 tubes/pack, 10 packs/case
801209	2.0 mL	Clear	Clear	100 tubes/pack, 10 packs/case
801275	0.5 mL	Amber	Amber	100 tubes/pack, 10 packs/case
801285	1.5 mL	Amber	Amber	100 tubes/pack, 10 packs/case
801295	2.0 mL	Amber	Amber	100 tubes/pack, 10 packs/case

### Screw Cap Micro Tubes Type A

P/N	Capacity	Tube Color	Cap Color	Configuration
801205A	0.5 mL	Clear	Clear	100 tubes/pack, 10 packs/case
801207A	1.5 mL	Clear	Clear	100 tubes/pack, 10 packs/case
801209A	2.0 mL	Clear	Clear	100 tubes/pack, 10 packs/case

### Optional colors for screw caps

Clear, Red, Blue, Green, Amber, Orange, Yellow, White, Violet, Pink